



MedicalDirector®  
**BLUECHIP**

## Day Surgery July 2017 Update for Australian Capital Territory MedicalDirector Blue Chip

Dear MedicalDirector Blue Chip Subscriber,

MedicalDirector is pleased to announce the release of the July 2017 update for MedicalDirector Blue Chip Day Surgery.

MedicalDirector strongly recommends that you apply updates when they become available – this update is valid from July 2017.

Please contact MedicalDirector Customer Care on **1300 300 161** (option 2) to schedule an installation time for this update.

### Included in this Update


---

- ICD-10-AM tables updated to Version 10
- Changes to HCP based extracts to reflect new ICD-10-AM Version 10
- Grouper interface upgraded for AR-DRG V9.0 and 3M CGS Version 2017.1
- Added new Health Fund myOwn to the list of health funds

## Important Notes on Installing this Update

### Installation Procedure:

1. Backup all DSW files.
2. Replace DSA Data.mde with version 7.00
3. Replace DSA.mdb with version 10.00
4. Start DSA and from **Advanced > Relink** re-link the Episodes and Data databases.
5. From **Advanced > Edit Datasets > Local Health Funds** add new entries for the new health fund code MYO, name "myOwn":



6. Add the new health fund into the Grouped Health Funds from **Advanced > Edit Datasets > Grouped Health Funds** as members of AHS:

Service	Fund	Delete
AHS	MBP	Medibank Private
AHS	MYO	myOwn
AHS	NHB	Navy Health Ltd

MedicalDirector Head Office  
 Level 5, 477 Pitt Street,  
 Haymarket NSW 2000, Australia  
 Tel: 1800 622 678  
 Fax: +61 2 9906 8910  
[www.medicaldirector.com](http://www.medicaldirector.com)  
 A.B.N. 76 068 458 515  
 Copyright© 2017