



## MedicalDirector Blue Chip

Dear MedicalDirector Blue Chip Subscriber,

MedicalDirector is pleased to announce the release of the July 2018 update for MedicalDirector Blue Chip Day Surgery.

MedicalDirector strongly recommends that you apply updates when they become available – this update is valid from July 2018.

Please contact MedicalDirector Customer Care on **1300 300 161** (option 2) to schedule an installation time for this update.

### Included in this Update

---

- Redefined third gender option as “Other”
- Adds option for “Residential Psychiatric Facility” for transfers and referrals.
- Updated table of Queensland Hospitals used for transfers.

### Important Notes on Installing this Update

#### Installation Procedure:

1. Make backups of everything.
2. Replace DSQ Data.mde with the new version (8.01).
3. Replace DSQ.mde with the new version (12.01).
4. Open DSQ and relink Episodes and Data as needed

