



MedicalDirector Bluechip Day Surgery July 2019 Update for South Australia

We are pleased to announce the release of the July 2019 update for MedicalDirector Bluechip Day Surgery update for South Australia.

We strongly recommend that you apply the following update as this update is valid from July 2019.

If you have any issues, contact MedicalDirector Customer Support on **1300 300 161**


Included in this Update

- Updated ICD10-AM tables to V11
- Updated HCP based extracts to reflect new ICD10 version
- Fixed problem with HCP extracts including morphology codes and/or leaving blanks in extract
- Updated National Private Patients Hospital Claim Form for 2019
- Updated Same Day Certificate for 2019 legislation
- Updated Grouper interface to support 3M CGS 2019.1 and DRG version 10
- Corrected errors in the ISAAC 1827 format extract
- Removed dropdown for year in ISAAC Control log

Installation Process

How to Backup Day Surgery Files


1. Download and extract the current SA Day Surgery file 2019-20 (DS_15.03_SA.zip) file to your local desktop
2. Open the DSS folder. *Referred to as DSS* (folder location can be found in Bluechip)
Setup > Practice > Links > Day Surgery



3. Create a backup folder within the DSS folder
4. Copy the DSS, DSS Data and DSS Episodes and ISAAC Control Log into the backup folder and rename with old e.g. DSS_old
5. Copy and replace the existing database DSS.mdb and the ISAAC Control Log from the SA Day Surgery file 2019-20 (DS_15.03_SA.zip) folder

Relinking Day Surgery module

1. Open Day Surgery (if prompted to create Episode select No)




2. When prompted to Find Episode select Cancel



3. Relink Episodes
Go to Add-ins > Advanced > Relink Databases > Episodes
4. Relink DS Data
Go to Add-ins > Advanced > Relink Databases > DS Data

5. Show links

Go to Add-ins > Advanced > Relink Databases > Show Links



6. Go to Add-ins > Advanced > Settings and update the DRG version to 9.0 or 10.0 if required

Property	Setting	Comments
DefAttendingDoctor		Default Attending Doctor
defCountry	1100	Default value for Country
DisCancer	N	Disable cancer notification
DRGVersion	9.0	Default DRG Version

7. Close Day Surgery